DISPOSABLE COAGULATION ACCESSORIES INSTRUCTIONS FOR USE AND CARE

ENGLISH (EN).....	2
ČESKY (CS)	4
DANSK (DA).....	6
DEUTSCH (DE)	8
Ελληνικά (EL)	10
ESPAÑOL (ES)	12
SUOMI (FI).....	14
FRANÇAIS (FR)	16
ITALIANO (IT)	18
LIETUVIŠKAI (LT).....	20
LATVISKI (LV).....	22
NEDERLANDS (NL).....	24
NORSK (NO)	26
POLSKI (PL).....	28
PORTUGUÊS (PT)	30
SVENSKA (SV)	32
TÜRKÇE (TR).....	34

ENGLISH (EN)

INSTRUCTIONS FOR USE AND CARE-DISPOSABLE COAGULATION ACCESSORIES

This product is disposable, packaged **STERILE** within a Tyvek pouch. Multiple pouches are placed in a secondary packaging level with a chipboard box. Multiple chipboard boxes are packed in a tertiary cardboard shipping case. Prior to the expiration date, contents are sterile if package has not been opened or damaged. For questions or additional information on our complete line of bipolar accessories, please contact Kirwan Surgical Products (KSP) at (781) 834-9500 or on the web at www.ksp.com.

INTENDED USE: Disposable coagulation handpieces are electrosurgical devices designed to be used in soft tissue procedures. Disposable cords are designed to connect an electrosurgical device to an electrosurgical generator. The 10-4601 Disposable Irrigation Bipolar Cord is designed to be compatible with the Kirwan AURA™ Irrigating Generator as well as the Malis™ CMC® Bipolar System and similar bipolar irrigating systems. The 10-4602 and 10-4602J Disposable Irrigation Bipolar Cord is designed to be compatible with Malis™ 1000 Irrigating Bipolar systems. The 10-4603 Disposable Irrigation Bipolar Cord is designed to be compatible with the Kirwan POLARIS™ Irrigation Module. Please refer to the labeling to determine the particular connectors and specific use of cords. Please refer to Kirwan Surgical Products LLC catalog or website (www.ksp.com) for descriptions of individual items.

TARGET USER/PATIENT GROUP: There is no target patient group for this device; the device is used at the surgeon's preference.

INDICATIONS FOR USE: Disposable coagulation handpieces are indicated for use when coagulation of soft tissue is required. Disposable cords are indicated for use when the connection of a coagulation handpiece to an electrosurgical generator is required. Disposable irrigation cords are indicated for use when the connection of a coagulation handpiece to an electrosurgical generator is required, as well as when the delivery of irrigation solution to the surgical site is required.

CONTRAINDICATIONS: None known.

CLINICAL BENEFITS: The intended clinical benefit of this device is the effective coagulation of soft tissue.

DEVICE CHARACTERISTICS: This device delivers energy from an electrosurgical generator to coagulate soft tissue.

CAUTION:

- Federal (USA) law restricts this device to sale by or on the order of a physician.
- Do not use this product if the sterile package is damaged or if the seal is broken.
- Users of these devices should be trained in the use of surgical instruments, high frequency electrosurgical devices, their accessories, and other related equipment.
- Any use of this device for tasks other than its intended purpose will usually result in a damaged or broken device.
- Irrigating cords are designed for use with collapsible irrigation solution containers only, and for use with bipolar electrosurgical units only.
- The non-booted connection of irrigation cords should be considered **"NOT PROTECTED"** against ingress of liquid with harmful effects (CODE IPX0).
- Frequent or excessive bending of finer (23-27) Gauge pencils may result in a damaged pencil tip.
- Because of the variability of output voltages and odes from generator to generator, **DO NOT USE** this accessory with generator setting exceeding the following voltages:
 - ◆ 1100Vp-p bipolar output for bipolar pencils and bipolar adapters.
 - ◆ 1200Vp-p bipolar output for bipolar forceps and bipolar cords.
 - ◆ 5000Vp-p monopolar output for monopolar forceps.
 - ◆ 7000Vp-p monopolar output for monopolar cords and monopolar adapters.

Refer to the appropriate electrosurgical generator manual for indications and instructions on voltage output characteristics to ensure that all safety precautions are followed. If no high frequency output is delivered to the accessory handpiece when the generator's activating switch is pressed, check the cord connection with the device and with the generator. If proper function is still not achieved and the accessory handpiece and generator function are confirmed as sound, replace the cord and refer the questionable cord to qualified personnel for further evaluation.

WARNINGS AND PRECAUTIONS:

1. Connect **BIPOLAR** accessories to the **BIPOLAR** receptacle only, and **MONOPOLAR** accessories to the **MONOPOLAR** receptacle. Improper connection of accessories may result in inadvertent accessory activation or other potentially hazardous conditions.
2. Proper disposal of devices and sharps possibly contaminated with blood, tissue, or other potentially infectious material present a biological risk and must be discarded in a closable, leak-proof, puncture-resistant receptacle that is adequately labeled (i.e., color coding or symbology) for easy identification as biohazard waste. Improper disposal may lead to infection of users or patients who come into contact with the device.

SETUP AND USE:

Forceps/Pencils	Cords	Irrigation Cords
<ol style="list-style-type: none"> 1. Attach the sterile forceps/pencils to the sterile cord ensuring that the forceps/pencil pins are fully seated in the cord receptacles. This condition ensures that the connection is splash proof. 2. The cords to the surgical electrodes should be positioned in such a way that contact with the patient or other leads is avoided. Temporarily unused ACTIVE ELECTRODES should be stored isolated from the patient. 3. Connect the sterile cord to the generator supply only while it is in the "OFF" or "STANDBY" mode. Failure to do so may result in injury or electrical shock to the patient or operating room personnel. 4. At the lowest power setting, test the forceps/pencils by pressing the generator's activating switch. If the generator fails to activate, check the forceps/pencils connection with the cord. 5. Power setting guidelines may vary due to differences in surgical techniques, patients, electrodes and surgical set-up. Start at the lowest power setting and increase as necessary to achieve the desired clinical effect. 	<ol style="list-style-type: none"> 1. Attach the sterile cord to the sterile instrument, ensuring that the contact pins are fully seated in the cord receptacles. 2. Connect the sterile cord to the generator supply only while it is in the "OFF" or "STANDBY" mode. Failure to do so may result in injury or electrical shock to the patient or operating room personnel. 3. At the lowest power setting, test the forceps/pencils by pressing the generator's activating switch. If the generator fails to activate, check the forceps/pencils connection with the cord. 4. At the lowest power setting, test the forceps/pencils by pressing the generator's activating switch. If the generator fails to activate, check the forceps/pencils connection with the cord. 5. Power setting guidelines may vary due to differences in surgical techniques, patients, electrodes and surgical set-up. Start at the lowest power setting and increase as necessary to achieve the desired clinical effect. 	<ol style="list-style-type: none"> 1. Close the roller clamp on the tubing. 2. Remove the cap from the spike and insert into the collapsible irrigation solution container. 3. Place the irrigation solution container onto an I.V. pole. 4. Insert the section of silicone tubing into the irrigation pump module if one is used refer to pump manufacturer's instruction. 5. Remove the luer cap and connect the male luer connector of the Irrigation Cord to the bipolar forceps. 6. Connect the Irrigation Cord plug to the bipolar forceps. 7. Insert the banana pins into the bipolar generator. 8. At the lowest power setting, test the forceps/pencils by pressing the generator's activating switch. If the generator fails to activate, check the forceps/pencils connection with the cord. 9. Power setting guidelines may vary due to differences in surgical techniques, patients, electrodes and surgical set-up. Start at the lowest power setting and increase as necessary to achieve the desired clinical effect.

RETURNS: For any inquiries, complaints, or serious incidents related to the use of this device, contact your Kirwan Surgical Products (KSP) distributor or KSP directly in the event a product has to be returned for evaluation. Call KSP customer service at (781) 834-9500 for a Return Merchandise Authorization (RMA). All serious incidents should also be reported to the competent authority for medical devices of your State. When returning products, follow these instructions:

- ◆ **DO NOT RETURN USED/CONTAMINATED DISPOSABLE DEVICES.** KSP shall not accept product it deems contaminated, and a health hazard to its employees.
- ◆ Ship product in a sturdy shipping box, with sufficient soft packaging materials to protect them.
- ◆ Secure with heavy tape, clearly identify the box as a return with the RMA# on the outside to expedite the process.
- ◆ Ship to:

AURA™ is a registered trademark of Kirwan Surgical Products LLC.

POLARIS™ is a trademark of Kirwan Surgical Products LLC.

LIBERTY™ is a trademark of Kirwan Surgical Products LLC.

Malis™ is a trademark of Synergetics IP, Inc.

CMC® is a registered trademark of Johnson & Johnson Professional, Inc.

Kirwan Surgical Products LLC
180 Enterprise Drive
Marshfield, MA 02050
(781) 834-9500

For Regulatory Issues Only
Medical Product Service GmbH
Borgasse 20
35619 Braunfels, GERMANY

Accumed
Viale Serfontana 10
6834 Morbio Inferiore
SWITZERLAND

Manufacturer

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

Catalog Number

Batch Code

Do Not Restertize

Use By Date

Caution

Sterile Barrier System –
Sterilized by Irradiation

Medical Device

Single Use –
Do Not Reuse

Do Not Use if
Package is Damaged

Authorized Representative in
the European Union

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CE 1639

**BACK COVER INTENTIONALLY
LEFT BLANK**

IFU-KDCA-5 Rev A (DOC-14662) Ver. 1

Approved By:

[\(CO-6571\) E6571](#)

Description

E6571. IFU-KDCA-5 REV A language headers for 2 languages were on the wrong page. Language headers were also moved to alternate which side they were on depending on each sheet's position in the book. Greek EC REP symbol was fixed.

Justification

The IFU's have not been printed yet. These issues were found during first proofing prior to the release of Rev A of the IFU. We will make the change to Rev A now without a revision change.

Assigned To:	Initiated By:	Priority:	Impact:
Tom Snow	Tom Snow	Low	Minor

Version History:

Author	Effective Date	CO#	Ver.	Status
Tom Snow	July 8, 2024 8:55 AM EDT	CO-6571	1	Published
Tom Snow	May 28, 2024 9:37 AM EDT	CO-6558	0	Superseded