

14-5000



14-5002



14-5004




14-5011



14-5137

Disposable Bipolar Pencil

ITEM #	TOTAL LENGTH	DESCRIPTION	GAUGE	WORKING LENGTH	FEATURES
14-5000	-	Non-Stick, Straight, Beveled	18	-	<b>NS</b>
14-5000N	-	Non-Stick, Straight, Black	18	-	<b>NS</b>
14-5002	-	Non-Stick, Curved	18	-	<b>NS</b>
14-5004	-	Straight	20	-	
14-5006	-	Tapered Tip, Sharp	20/23	-	
14-5006N	-	Tapered Tip, Sharp, Black	20/23	-	
14-5007	-	Straight	23	-	
14-5007L	-	Straight, Long	23	-	
14-5009	-	Tapered Tip, Blunt	20/23	-	
14-5010	-	Widestroke Dual Coag Non-Stick, 30° Angle	18	-	<b>NS</b>
14-5011	-	Straight	25	-	
14-5011L	-	Straight, Long	25	-	
14-5027	-	Straight	27	-	
14-5137	-	Bent, 3mm Extended Tip	23	-	



14-7012



14-7013



14-7014



14-7016



14-7018



### Reusable TPR Bipolar Pencil

ITEM #	TOTAL LENGTH	DESCRIPTION	GAUGE	WORKING LENGTH	FEATURES
14-7000	-	Non-Stick, Straight	18	-	<b>NS</b>
14-7002	-	Non-Stick, Curved	18	-	<b>NS</b>
14-7006	-	Straight	20	-	-
14-7007	-	Straight	23	-	-
14-7008	-	Tapered Tip, Sharp	20/23	-	-
14-7009	-	Tapered Tip, Blunt	20/23	-	-
14-7010	-	Widestroke Dual Coag Non-Stick, 30° Angle	18	-	<b>NS</b>
14-7012	7" (17.8cm)	Sinusoidal Non-Stick, Straight	14	3" (7.6cm)	<b>NS</b>
14-7013	-	10mm Extended Center Tip	20	-	-
14-7014	7" (17.8cm)	Sinusoidal Non-Stick, Angled	14	3" (7.6cm)	<b>NS</b>
14-7016	8" (20.3cm)	Non-Stick, 10° Angled, Pencil Point Tip	14	4" (10.2cm)	<b>NS</b>
14-7017	12" (30.5cm)	Non-Stick, Dome Tip	14	8" (20.3cm)	<b>NS</b>
14-7018	15 3/4" (40cm)	Non-Stick, 10° Angled, Pencil Point Tip	14	11 3/4" (29.8cm)	<b>NS</b>