




## DISPOSABLE FLYING LEAD BIPOLAR CORDS

ITEM #	LENGTH	DESCRIPTION	FEATURES
10-4000	12' (3.6m)	Disposable bipolar cord sterile green	②
10-4001	12' (3.6m)	Disposable bipolar cord sterile clear	②
10-4001A	12' (3.6m)	Disposable bipolar cord sterile clear, 90 degree angle	②
10-4021	16.4' (5m)	Disposable bipolar cord sterile clear, bladed connector, sterile clear	②
10-4601	12' (3.6m)	Disposable bipolar irrigation tube and cord assembly, sterile clear	②
10-4602	12' (3.6m)	Disposable bipolar irrigation tube and cord assembly, M1000 type, sterile clear	②
10-4611	12' (3.6m)	Disposable bipolar dual irrigation cord assembly, sterile clear	②
10-4603	12' (3.6m)	Disposable bipolar irrigation tube and cord assembly, POLARIS, sterile clear	② PL
10-4601A	12' (3.6m)	Disposable bipolar irrigation tube and cord assembly, 90 degree angle, sterile clear	②
10-4602A	12' (3.6m)	Disposable bipolar irrigation tube and cord assembly, M1000 type, 90 degree angle, sterile clear	②



## REUSABLE FLYING LEAD BIPOLAR CORDS

ITEM #	LENGTH	DESCRIPTION
10-5000	12' (3.6m)	Reusable TPR bipolar cord (thermoplastic rubber)
10-5000A	12' (3.6m)	Reusable TPR bipolar cord (thermoplastic rubber), 90 degree angle
10-5001	6' (1.8m)	Reusable TPR bipolar cord (thermoplastic rubber)
10-6000	12' (3.6m)	Reusable silicone bipolar cord




## REUSABLE FLYING LEAD BIPOLAR CORDS

ITEM #	LENGTH	DESCRIPTION
10-6002	12' (3.6m)	Reusable silicone bipolar scissors cord



## REUSABLE BIPOLAR CORDS

ITEM #	LENGTH	DESCRIPTION
10-6010	12' (3.6m)	Reusable silicone bipolar cord with Lemo plug (for Storz-type Daisy)



10-6011	12' (3.6m)	Reusable silicone bipolar cord with Lemo plug (for Storz-type Premier)
10-6012	12' (3.6m)	Reusable silicone bipolar cord with 2 Lemo plugs



10-6020	12' (3.6m)	Reusable silicone bipolar cord (with banana pins & bladed connector)
---------	------------	--



10-6025	12' (3.6m)	Reusable silicone bipolar cord with booted forceps connector (for Wolf/Martin-type)
---------	------------	---



10-6037	12' (3.6m)	Reusable silicone bipolar cord with Erbe-type, EU connector
---------	------------	---



# REUSABLE BIPOLAR CORDS

ITEM #      LENGTH      DESCRIPTION

10-6050    12' (3.6m)    Reusable silicone bipolar cord with Erbe-type/Storz Millennium-type booted forceps connector



10-6055    12' (3.6m)    Reusable silicone bipolar cord with Redel connector



10-6060    12' (3.6m)    Reusable silicone bipolar cord with Fischer plug



10-6075    12' (3.6m)    Reusable silicone bipolar cord (for O.M.S. -type Phaco machine)

10-6076    6' (1.8m)    Reusable silicone bipolar cord (for O.M.S. -type Phaco machine)















10-7001    12' (3.6m)    Reusable silicone bipolar handswitch cord



10-7050    6'4" (1.9m)    Reusable blue silicone coaxial cord (for D.O.R.C. -type Phaco)




# DISPOSABLE FIXED PIN BIPOLAR CORDS

ITEM #	LENGTH	DESCRIPTION	FEATURES	
10-4001J	12' (3.6m)	Disposable bipolar cord with fixed pins, 38mm, sterile clear		
10-4001R	12' (3.6m)	Disposable bipolar cord with fixed pins, 22mm, sterile clear		
10-4000V	12' (3.6m)	Disposable bipolar cord with fixed pins, 28.6mm, sterile green		
10-4001V	12' (3.6m)	Disposable bipolar cord with fixed pins, 28.6mm, sterile clear		
10-4601J	12' (3.6m)	Disposable bipolar irrigation tube and cord assembly with fixed pins, 38mm, sterile clear		
10-4602J	12' (3.6m)	Disposable bipolar irrigation tube and cord assembly, M1000 type with fixed pins, 38mm, sterile clear		
10-4601V	12' (3.6m)	Disposable bipolar irrigation tube and cord assembly with fixed pins, 28.6mm, sterile clear		

# REUSABLE FIXED PIN BIPOLAR CORDS

ITEM #	LENGTH	DESCRIPTION
10-5000R	12' (3.6m)	Reusable TPR bipolar cord with fixed pins, 22mm
10-6000R	12' (3.6m)	Reusable silicone bipolar cord with fixed pins, 22mm
10-5000V	12' (3.6m)	Reusable TPR bipolar cord with fixed pins, 28.6mm
10-6000V	12' (3.6m)	Reusable silicone bipolar cord with fixed pins, 28.6mm
10-6000J	12' (3.6m)	Reusable silicone bipolar cord with fixed pins, 38mm
10-6009	12' (3.6m)	Reusable silicone bipolar cord with flexible molded plug
10-6020J	12' (3.6m)	Reusable silicone bipolar cord with fixed pins, 38mm, EU connector
10-6029	12' (3.6m)	Reusable silicone bipolar cord with flexible molded plug, EU connector



ITEM # DESCRIPTION

10-6500 Reusable bipolar Erbe-type/Storz-Millennium-type adapter (Erbe-type to banana pin)



10-6575 Reusable O.M.S.-type Phaco machine adapter



10-6580 Reusable D.O.R.C. type Phaco machine adapter



10-6585 Reusable Geuder Megatron machine adapter



10-6590 Reusable Oertli adapter



10-6595 Reusable Ellman-type adapter 4.0mm to 3.0mm banana pin



30-5110 Reusable monopolar suction coagulator adapter (banana pin to Bovie)

30-5205 Reusable monopolar suction coagulator adapter (resectoscope to Bovie)



30-5210 Reusable monopolar suction coagulator adapter (resectoscope to banana pin)



# DISPOSABLE MONOPOLAR CORDS

ITEM #	LENGTH	DESCRIPTION	FEATURES
30-1201	10' (3m)	Disposable monopolar endoscopic cord, sterile green with straight plug	
30-1201L	10' (3m)	Disposable monopolar endoscopic cord, sterile green with extended straight plug	
30-1203	10' (3m)	Disposable monopolar endoscopic cord, sterile green with 90 degree angled plug	
30-1204	10' (3m)	Disposable monopolar endoscopic cord, sterile green with banana pin	
30-1208	10' (3m)	Disposable monopolar endoscopic cord, sterile green with Luer connector	
30-1209	10' (3m)	Disposable monopolar endoscopic cord, sterile green with 2mm pin	
30-1211	10' (3m)	Disposable monopolar endoscopic cord, sterile green with 3mm pin	



# REUSABLE MONOPOLAR CORDS

ITEM #      LENGTH      DESCRIPTION

30-1000    10' (3m)    Reusable silicone monopolar cord with Kirwan plug



30-1205    10' (3m)    Reusable silicone monopolar endoscopic cord with Kirwan plug



30-1206    10' (3m)    Reusable silicone monopolar endoscopic cord with Kirwan plug and banana pin



30-3003    10' (3m)    Reusable silicone monopolar cord with banana pins





# INTERNATIONAL BIPOLAR CORDS

ITEM #      DESCRIPTION      FEATURES

10-4018J Disposable 5m bipolar cord, sterile, fixed 38mm pins to European bladed connector



10-4021V Disposable 5m bipolar cord, sterile, fixed 28.6mm pins to European bladed connector



10-4028J Disposable 5m bipolar cord, sterile, fixed 38mm pins to booted dual connector




10-4038 Disposable 5m bipolar cord, sterile, Erbe-type to European bladed connector



10-4048 Disposable 5m bipolar cord, sterile, Wolf/Martin-type to European bladed connector



10-4058 Disposable 5m bipolar cord, sterile, Erbe-type to booted dual connector



10-4068 Disposable 5m bipolar cord, sterile, Wolf/Martin-type to booted dual connector

