




Disposable Tungsten Microdissection Needle, Insulated, Straight

| ITEM # | TOTAL LENGTH | DESCRIPTION | FEATURES |
|---------|-----------------|-------------|----------|
| 30-8000 | 1 1/4" (3.2cm) | - | IN |
| 30-8002 | 2" (5.1cm) | Extra Fine | IN |
| 30-8004 | 5 1/2" (14cm) | - | IN |
| 30-8006 | 1 5/8" (4.1cm) | - | IN |
| 30-8008 | 2" (5.1cm) | - | IN |
| 30-8010 | 2 1/2" (6.4cm) | - | IN |
| 30-8011 | 4 5/8" (11.7cm) | - | IN |
| 30-8013 | 6 5/8" (16.8cm) | - | IN |



| ITEM # | TOTAL LENGTH | DESCRIPTION | FEATURES |
|---------|--------------|-------------------|----------|
| 30-8015 | 4" (10.2cm) | Insulated, Curved | IN |



Disposable Tungsten Microdissection Needle, Insulated, Angled 45 Degrees

| ITEM # | TOTAL LENGTH | DESCRIPTION | FEATURES |
|---------|-----------------|-------------|----------|
| 30-8007 | 1 5/8" (4.1cm) | - | IN |
| 30-8009 | 2" (5.1cm) | - | IN |
| 30-8012 | 4 1/2" (11.4cm) | - | IN |
| 30-8014 | 6 1/2" (16.5cm) | - | IN |
| 30-8017 | 4 7/8" (12.4cm) | - | IN |
| 30-8018 | 6 7/8" (17.5cm) | - | IN |



| ITEM # | TOTAL LENGTH | DESCRIPTION | FEATURES |
|---------|-----------------|-------------------|----------|
| 30-8016 | 5 7/8" (14.9cm) | Insulated, Curved | IN |

30-6011



Disposable Monopolar Electrode, Insulated

| ITEM # | TOTAL LENGTH | TIP ROD DIAMETER | INSULATION DIAMETER | DESCRIPTION | FEATURES |
|---------|--------------|------------------|---------------------|-------------|----------|
| 30-6011 | 30.2cm | 1.37mm | 1.72mm | Tapered Tip | IN |
| 30-6012 | 30.2cm | 1.37mm | 1.72mm | Blunt Tip | IN |
| 30-6020 | 25cm | 1.37mm | 1.72mm | Blunt Tip | IN |
| 30-6021 | 25cm | 1.32mm | 1.67mm | Tapered Tip | IN |
| 30-6022 | 25cm | 1.32mm | 1.67mm | Dome Tip | IN |
| 30-6023 | 25cm | 1.00mm | 1.32mm | Dome Tip | IN |

Electrode Tip Detail



Disposable Monopolar Suction Stimulator

| ITEM # | WORKING LENGTH | DESCRIPTION | FEATURES |
|---------|----------------|-------------|----------|
| 41-9012 | 8cm | 5 French | IN |
| 41-9014 | 10cm | 5 French | IN |
| 41-9016 | 13cm | 5 French | IN |
| 41-9022 | 8cm | 7 French | IN |
| 41-9024 | 10cm | 7 French | IN |
| 41-9026 | 13cm | 7 French | IN |
| 41-9032 | 8cm | 10 French | IN |
| 41-9034 | 10cm | 10 French | IN |
| 41-9036 | 13cm | 10 French | IN |
| 41-9042 | 8cm | 12 French | IN |
| 41-9044 | 10cm | 12 French | IN |
| 41-9046 | 13cm | 12 French | IN |

30-6001



Reusable Monopolar Electrode, Insulated

| ITEM # | TOTAL LENGTH | DESCRIPTION | FEATURES |
|---------|--------------|-------------|----------|
| 30-6001 | 15" (38.1cm) | Hook Tip | IN |
| 30-6002 | 15" (38.1cm) | Spoon Tip | IN |
| 30-6003 | 15" (38.1cm) | Spatula Tip | IN |
| 30-6004 | 15" (38.1cm) | Knife Tip | IN |

30-4003



Reusable Monopolar Chuck Handle

| ITEM # | TOTAL LENGTH | DESCRIPTION | FEATURES |
|---------|-----------------|---------------|----------|
| 30-4000 | 6" (15.2cm) | With 6' Cord | - |
| 30-4001 | 6" (15.2cm) | With 10' Cord | - |
| 30-4002 | 6" (15.2cm) | With 15' Cord | - |
| 30-4003 | 8 1/2" (21.6cm) | With 6' Cord | - |
| 30-4004 | 8 1/2" (21.6cm) | With 10' Cord | - |
| 30-4005 | 8 1/2" (21.6cm) | With 15' Cord | - |

30-4100



Reusable Monopolar Cut & Coag Push Button Pencil

| ITEM # | TOTAL LENGTH | DESCRIPTION | FEATURES |
|---------|--------------|---------------|----------|
| 30-4100 | - | With 10' Cord | - |